

M95

## Simulator Pediatric IV Hand



### Set of pediatric hands for practicing IV

Become confident in safety insertion of intravenous cannulas on small infant hands with thin and soft veins

#### FEATURES

- 1 | Fingers and dorsum of hand are flexible, able to be gripped, and
- 2 | repositionable as in real-life procedures.
- 3 | Two puncture sites on the back of the hand
- 4 | Practice dressing and securing the cannula

#### SKILLS

- | Pediatric dorsal veins injection
- | Blood collection
- | Intravenous infusion

#### DESCRIPTIONS

##### SPECIFICATIONS

Size: W23 x D6 x H5 cm/ W9 x D2.3 x 1.9 inch (one-year-old) Power: AC100-240V  
W30 x D8 x H6 cm/ W11.8 x D3.1 x 2.3 inch (three-year-old) Power Consumption: 5W

##### SET INCLUDES

2 injection arm-hands (one-year-old)	1 spoon
2 injection arm-hands (three-year-old)	1 plastic jar
2 stands	1 plastic beaker
1 circulation pump	1 storage case
1 jar of coloring powder (red)	1 instruction manual

##### MATERIALS

Soft resin, hard resin  
Latex free

##### REPLACEMENT PARTS

11350-010	1 replacement injection arm-hand (three-year-old)
11350-020	1 replacement injection arm-hand (one-year-old)

##### RECOMMENDED DEVICES

indwelling needle: 23G  
butterfly needle: 23G

\* Needle not included.

M43D


## Pediatric Lumbar Puncture Simulator II

### Acquire skills for safe pediatric LP

To keep the patient body safely at the appropriate position is critical in making successful lumbar puncture. This full-body soft and bendable baby manikin facilitates training with reality.



Level of posterior superior iliac crests



#### FEATURES

- 1 | Soft, seamless baby manikin allows practice in safe patient
- 2 | positioning and holding.
- 3 | Two types of puncture pads: deep or shallow
- 4 | Realistic needle tip resistance
- 5 | Anatomically correct landmarks for palpation

#### SKILLS

- | Patient positioning
- | Palpation of landmarks
- | Lumbar puncture and CSF collection
- | CSF pressure measurement

#### Product Supervision

Takahiro Amano, M.D.  
Vice President  
Senior Vice-Dean, Postgraduate School  
Professor and Director, Center of Postgraduate  
Medical Education International University of Health  
and Welfare,  
Honorary Director, Sanno Medical Center  
Gregory A. Plotnikoff, M.D., M.T.S., F.A.C.P.  
Medical Director, Abbott Northwestern Hospital

#### DESCRIPTIONS

##### SPECIFICATIONS

Size: W43 x D34 x H20 cm Weight: 1.5 kg/ 3.3 lbs  
W17 x D13.4 x H7.8 inch

##### SET INCLUDES

1 baby manikin (7-10 months)	1 syringe (50mL)
2 pediatric LP block (shallow)	1 instruction manual
2 pediatric LP block (deep)	
1 irrigation bag	
1 stand for irrigation bag	

##### MATERIALS

Soft resin  
Latex free

##### REPLACEMENT PARTS

11348-100-04	stand for irrigation bag
11348-100-05	irrigation bag
11348-010	pediatric LP block (shallow)
11348-420	pediatric LP block (deep)

##### RECOMMENDED DEVICES

22G-24G spinal needle

\* Needles are not included.