

MW1



# Pitting Edema Simulator

## Five variations of pitting edema simulator including a normal one

Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"

### Product Supervision

Junji Otaki  
Professor of Tokyo Medical University



### FEATURES

- 1 | Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"
- 2 | The models can be attached to the lower leg of a simulated patients or a training manikins.
- 3 | 4 levels of pitting edema according to pressures and the depth of indentations.
- 4 | Repeatable and accurate, as if palpating actual skin surface to the fibula.

### SKILLS

- | Differentiation of the pitting edema: four different levels of edemas and a normal case
- | Assessment of edemas on the lower legs

### KEY FEATURES



Normal

1+ Light Pitting

2+ Moderate Pitting

3+ Deep Pitting

4+ Very Deep Pitting

### DESCRIPTIONS

#### SPECIFICATIONS

Size: 10x23cm/3.9x9 inch

#### SET INCLUDES

5 variations of strap-on pitting edema model 1 piece each

#### MATERIALS

Soft resin  
Latex free